

2013.10.23

**Risø HDRL Technical note E3****Evaluation of dose map for irradiation at verification dose  $D_{\text{ver}}$ .**

Author: Arne Miller  
Risø High Dose Reference Laboratory  
DTU Nutech  
DK-4000 Roskilde  
Denmark

The standard for radiation sterilization of medical devices EN / ISO 11137, part 2 specifies how product samples shall be irradiated at verification dose,  $D_{\text{ver}}$ , in connection with establishing, substantiating or auditing a dose for sterilization,  $D_{\text{ste}}$ . The procedure used for sterilization dose establishment for Method 1 as described in the standard is used as an example in this note. The same principles can be applied for Method 2, and Method  $VD_{\text{max}}$ .

**EN / ISO 11137, part 2:****7.2.5.2**

Irradiate the product items at the verification dose. Determine the dose. If the highest dose to the product items exceeds the verification dose by more than 10 %, and the sterilization dose is to be established using Method 1, the verification dose experiment shall be repeated. If the arithmetic mean of the highest and lowest doses delivered to product items is less than 90 % of the verification dose, the verification dose experiment may be repeated. If this mean dose is less than 90 % of the verification dose and, on performance of the tests of sterility, acceptable results are observed (see 7.2.6.1), the verification experiment need not be repeated.

In case of gamma irradiation it is normal to measure the delivered doses during irradiation of product samples, and obtain maximum and minimum doses to product samples. In case of e-beam irradiation dose mapping can be carried out in order to determine a monitoring dose  $D_{\text{mon}}$  that shall be used for irradiation of product samples in order irradiate these at  $D_{\text{ver}}$ . Measurement uncertainties must be taken into account in order to determine the required  $D_{\text{mon}}$  and to fulfill the requirement above.

Electron accelerator irradiation of product samples is usually carried out in a geometry that is designed to minimize the difference between maximum and minimum dose to product. The geometry is often selected to be single layer of product samples placed between two 3-mm polystyrene plates. This package is irradiated from two sides. The polystyrene plates serve to scatter the e-beam leading to a more uniform dose distribution and also to limit backscatter from metal components of conveyor or carrier.

Dose mapping is usually carried out by irradiating with doses that are conveniently measured, and based on the results a required monitoring dose  $D_{\text{mon}}$  for irradiation of product samples at  $D_{\text{ver}}$  can be calculated.

For e-beam irradiation it can be expected that there is significant difference between maximum and minimum doses, even within a single product sample. It might therefore be needed to place several dosimeters at the product samples in order to obtain valid dose map results. The results of the dose map can be expressed in terms of average minimum dose  $D_{\text{min}}$  and average maximum dose  $D_{\text{max}}$  at each of the individual dosimeter positions accompanied by their respective uncertainties. The minimum dose  $D_{\text{min}}$  and the maximum dose  $D_{\text{max}}$ , which shall be considered with respect to the requirements of the standard, are the lowest and the highest of the individual average minimum and maximum doses, respectively.

The mean dose is calculated as  $D_{\text{mean}} = (D_{\text{max}} + D_{\text{min}}) / 2$ .

The required  $D_{\text{mon}}$  should be determined to be in a range that will allow  $D_{\text{mean}}$  to be more than 90% of  $D_{\text{ver}}$ , while  $D_{\text{max}}$  does not exceed  $D_{\text{ver}}$  by more than 10%.

In order not to overdose product samples, the upper 95% confidence limit of  $D_{\text{max}}$  (DU) is used to calculate the maximum limit of  $D_{\text{mon}}$ .

The maximum limit of  $D_{\text{mon}}$  is calculated as

$$D_{(\text{mon}, \text{max})} = D_{(\text{mon}, \text{meas})} * D_{\text{ver}} * 1.1 / \text{DU}$$

The minimum limit of  $D_{\text{mon}}$  is calculated as

$$D_{(\text{mon}, \text{min})} = D_{(\text{mon}, \text{meas})} * D_{\text{ver}} * 0.9 / D_{\text{mean}}$$

An Excel file is attached that shows an example for calculations of the required range of monitor dose  $D_{\text{mon}}$ .

RISØ HIGH DOSE REFERENCE LABORATORY  
 Template 05 (1) Risø B3 / RisøScan

Date: 2013.10.17

Report number: 13R-xx  
 Customer: I.D.  
 Product: product i.d.  
 Facility: E-beam, 2-sided  
 Irradiation for Dver  
 Irradiated between PS plates

Operator: tk  
 Env. Temp: 22  
 Env. R.H.: 37%

RisøScan Data

Alan. Ref.AR602/1140	14.84	0.993 Correction factor
B3 ref	14.95	

BOX 1

	D(min)		D(max)	
	Uncorrected	Corrected	Uncorrected	Corrected
1.1	25.5	25.3	28.2	28.0
1.2	27.4	27.2	29.8	29.6
1.3	26.8	26.6	29.5	29.3
1.4	26.1	25.9	30.0	29.8
2.1	25.3	25.1	28.1	27.9
2.2	26.2	26.0	29.2	29.0
2.3	25.8	25.6	29.5	29.3
2.4	25.4	25.2	29.9	29.7
3.1	24.5	24.3	28.1	27.9
3.2	25.8	25.6	29.4	29.2
3.3	26.3	26.1	30.2	30.0
3.4	25.8	25.6	30.4	30.2
4.1	25.7	25.5	27.7	27.5
4.2	25.6	25.4	29.3	29.1
4.3	26.7	26.5	29.9	29.7
4.4	26.0	25.8	30.9	30.7
5.1	25.7	25.5	28.6	28.4
5.2	26.6	26.4	29.0	28.8
5.3	26.4	26.2	29.8	29.6
5.4	26.6	26.4	29.6	29.4

B3 ref measurement	
1	15
2	14.9
Average	14.95

Dmon, requested 2 x	15	kGy
Dmon, measured 2x	14.8	kGy

Average position #	D(min)			D(max)		
	Average kGy	s.d. kGy	%	Average kGy	s.d. kGy	%
1	25.2	0.5	2.0	27.9	0.3	1.1
2	26.1	0.7	2.7	29.1	0.3	1.0
3	26.2	0.4	1.5	29.6	0.3	1.0
4	25.8	0.4	1.7	29.9	0.5	1.7

Dmin	25.2	0.5	2.0	Dmax	29.9	0.5	1.7
D(mean)	27.5			D(max)/Dmin	1.19		

Uncertainties (k=1), %			
	D(min)	D(max)	Dmon
Dose map	2.0	1.7	2.3
Calibration	2.1	2.1	
Facility var.	1	1	
at k = 1, %	3.8	3.7	
at k = 1, kGy	0.96	1.10	

Upper 95% conf. Limit of D(max), DU	32.1	kGy
Dver, kGy	7.0	Insert value based on bioburden data

Dmon(max), 2x	3.55	kGy	Dmon(max) = D(mon, meas)*Dver*1.1/DU
Dmon(min) 2x	3.38	kGy	Dmon(min) = D(mon, meas)*Dver*0.9/Dmean
Required Dmon = (Dmon(max)+Dmon(min)) / 2 =	3.47	kGy	

The maximum dose to product may be 10% greater than the verification dose  
 The mean dose to product may be 90% of the verification dose